



**Feel Well**  
MEDICAL



# PRODUCT CATALOGUE



[www.feelwellmedical.com](http://www.feelwellmedical.com)



**Feel Well**  
MEDICAL

## Bilateral Inguinal Truss

In cases of inguinal hernia, it can be easily used under underwear. Its application before surgery prevents the growth of the hernia. Used to support the surgical area after surgery. Thanks to its special scrotal pad, when worn it does not disturb the body during use. Has arches that allow it to align and maintain its proper position.



FW-100

## Postpartum Corset

Supports the postpartum belly with its designed elastic tissue. Prevents skin and organ sagging. Supports firming and supports upright posture, especially during breastfeeding.



FW-101

## Elastic Abdominal Corset

Made from premium fabric that does not sweat and is breathable. Can be used easily without being visible under the clothes. Supports the belly and surgical area with its designed elastic texture."



FW-102

## Elastic Abdominal Corset (Standard)

Made from premium fabric that does not sweat and is breathable. Can be used easily without being visible under the clothes. Supports the belly and surgical area with its designed elastic texture."



FW-103

## Umbilical Hernia Corset

Includes silicone pad covered with breathable cotton. With its designed elastic tissue, it supports the hernia area and prevents its growth.

FW-104



## Umbilical Hernia Corset (Standard)

Includes silicone pad covered with breathable cotton. With its designed elastic tissue, it supports the hernia area and prevents its growth.

FW-105



## Colostomy Corset (16cm)

Reduces the risk of postoperative parastomal hernia. Supports stoma after surgery. Protect the colostomy bag. Made of premium fabric that does not sweat and is breathable. Can be used easily without being visible under clothing. With two removable plastic underwire supports on both sides of the stoma opening, it prevents unwanted movements of the stoma and possible leaks that may occur during daily sports activities.

FW-106



## Colostomy Corset (26cm)

Reduces the risk of postoperative parastomal hernia. Supports stoma after surgery. Protect the colostomy bag. Made of premium fabric that does not sweat and is breathable. Can be used easily without being visible under clothing. With two removable plastic underwire supports on both sides of the stoma opening, it prevents unwanted movements of the stoma and possible leaks that may occur during daily sports activities.

FW-107



## Neoprene Abdominal Corset

Made from neoprene fabric. Supports the navel and surgical area. By helping to remove abdominal and belly fat, it provides slimming and tightening in the waist area of the users. With regular use and a healthy eating routine, it allows people to have the desired appearance and body lines.



FW-108

## Unilateral Inguinal Truss

In cases of inguinal hernia, it can be easily used under underwear. Its application before surgery prevents the growth of the hernia. Used to support the surgical area after surgery. Thanks to its special scrotal pad, when worn it does not disturb the body during use. Has arches that allow it to align and maintain its proper position.



FW-109

## Working Corset

It is 16 cm high. Made of fabric that is breathable and supports the waist and back area. There are 4 plastic supports in the back area. Stabilization is provided with two cross-band belts from the back area to the waist area.

FW-110



## Wool Camel Hair Corset

Made from camel hair fabric. Provides warmth and compression to the abdomen and waist area. It provides warmth and compression by preserving body temperature in the applied organ.

FW-111



## Pregnancy Corset

It is made from breathable cotton and soft fabric. It can be used before and after birth. It prevents waist and back pain by supporting the abdomen during pregnancy. It is also used as an abdominal tightener after birth. It is compatible with the growth of the abdomen with its adjustable straps.

FW-112



## Lumbostad Corset With Steel Underwire (32 Cm)

Made from flexible, breathable, anti-sweat fabric and special elastic material. The 32cm height helps to support the lumbosacral region. Consists of two steel and four plastic underwire supports that stabilize the spine. The steel underwire is removable and removable. The steel underwires can be bent manually according to lumbar lordosis. Reduces the load on the discs with its cross-elastic band belt. Elastic fabric provides compression to the waist area.

FW-113



## Lumbostad Corset With Steel Underwire (32 Cm) (Standard)

Made from flexible, breathable, anti-sweat fabric and special elastic material. The 32cm height helps to support the lumbosacral region. Consists of two steel and four plastic underwire supports that stabilize the spine. The steel underwire is removable and removable. The steel underwires can be bent manually according to lumbar lordosis. Reduces the load on the discs with its cross-elastic band belt. Elastic fabric provides compression to the waist area.

FW-114



## Lumbostad Corset With Steel Underwire (26 Cm)

Made of flexible, breathable and anti-sweat fabric. The height is 26cm. Consists of two steel and four plastic underwire supports that stabilize the spine. The steel underwire is removable and attachable. Steel underwires can be bent manually according to lumbar lordosis feature. Reduces the load on the discs with its cross-elastic band belt. Elastic fabric provides compression to the waist area.

FW-115



## Lumbostad Corset With Steel Underwire (26 Cm) (Standard)

Made of flexible, breathable and anti-sweat fabric. The height is 26cm. Consists of two steel and four plastic underwire supports that stabilize the spine. The steel underwire is removable and attachable. Steel underwires can be bent manually according to lumbar lordosis feature. Reduces the load on the discs with its cross-elastic band belt. Elastic fabric provides compression to the waist area.

FW-116



## Woolen Corset

Made from elastic wool fabric. Provides warmth and compression to the abdomen and waist area. It provides warmth and compression by preserving body temperature in the applied organ.

FW-117



## Dorso-Lumbar Corset

Made from lightweight, breathable, flexible, moisture-absorbing cotton blended elastic fabric. The structure of the fabric and the design of the corset ensure that the pressure is distributed evenly. The design of the support promotes proper posture by pulling the shoulders back and aligning the spine. Consists of two steel and four plastic underwire supports. Removable metal supports provide additional stable support. The soft pad at the waist keeps it tension-free as well as providing extra stability. Adjustable cross straps ensure individual fit and support.



FW-118



## Dorso-Lumbar Corset (Standard)

Made from lightweight, breathable, flexible, moisture-absorbing cotton blended elastic fabric. The structure of the fabric and the design of the corset ensure that the pressure is distributed evenly. The design of the support promotes proper posture by pulling the shoulders back and aligning the spine. Consists of two steel and four plastic underwire supports. Removable metal supports provide additional stable support. The soft pad at the waist keeps it tension-free as well as providing extra stability. Adjustable cross straps ensure individual fit and support.



FW-119



## Postorex With Steel Underwire Corset (Standard)

Made from breathable and sweatproof fabric. There are steel underwires on both sides of the spine. Prevents hyperkyphotic dysphagia and its progression.

FW-120



## Steel Supported Chest Corset

Has sternum support with steel splints and is made of high-quality elastic fabric. It is easy to put on and take off thanks to the practical Velcro system and adjustable shoulder straps. Supports the sternum and thoracic cranium and stabilizes the sternoclavicular joints and rib cage after open heart and thoracic surgeries. Provides protection against pain caused by sudden movements and gives the patient confidence and comfort, Helps recovery after surgery with its proprioceptive effect.

FW-121



## Hyperextension Corset

Its outer surface is made of special aluminum material. The part in contact with the skin surface is covered with a terry pad. It is a body orthosis that stimulates hyperextension position by limiting the flexion of the trunk from 3 points. Length and width adjustments can be made on the orthosis thanks to the screw system on it. Provides anatomically correct posture by limiting spinal movements to a certain extent. Accelerates the healing process by increasing stability.

FW-122



## Posturex Bandage

Made from breathable and sweatproof fabric. There are plastic underwires on both sides of the spine. Prevents hyperkyphotic dysphagia and its progression.

FW-123



## Postorex Bandage (Standard)

Made from breathable and sweatproof fabric. There are plastic underwires on both sides of the spine. Prevents hyperkyphotic dysphagia and its progression.

FW-124



## Postorex Wide Underwire Corset (Standard)

Made from neoprene fabric. There are plastic underwires extending to the lumbar region on both sides of the spine. Belts at the waist position provide support for the lumbar cavity. It prevents hyperkyphotic posture disorder and has a stimulating effect.

FW-125



## Postorex Single Underwire (Standard)

Made from neoprene fabric. Has an underwire that will support the spine from the middle. Prevents Hyper Kyphotic posture disorder and provides a stimulating effect.

FW-126



## Postorex Long Underwire Corset (Standard)

Made from neoprene fabric. There are plastic underwires extending to the lumbar region on both sides of the spine. Belts at the waist position provide support for the lumbar cavity. It prevents hyperkyphotic posture disorder and has a stimulating effect.

FW-127



## Calf Support

Made from neoprene material. Increases blood flow, supports the calf and increases temperature. Reduces edema with compression and provides an analgesic effect. It is suitable for use in calf muscle injuries.



FW-128

## Calf Support (Standard)

Made from neoprene material. Increases blood flow, supports the calf and increases temperature. Reduces edema with compression and provides an analgesic effect. It is suitable for use in calf muscle injuries.



FW-129

## Thigh Support

Made from neoprene material. Provides analgesic effect with compression. Reduces edema and supports the muscles. Used for injuries of abductor, adductor, hamstring and quadriceps muscles.



FW-130

## Thigh Support (Standard)

Made from neoprene material. Provides analgesic effect with compression. Reduces edema and supports the muscles. Used for injuries of abductor, adductor, hamstring and quadriceps muscles.



FW-131

## Sling With Abduction Pillow (30 Degree)

It is easy to use and practical thanks to its comfortable and adjustable straps. Designed with a pad to keep the shoulder at 30-45 degrees.

FW-132



## Sling With Abduction Pillow (45 Degree)

It is easy to use and practical thanks to its comfortable and adjustable straps. Has a pad that keeps the shoulder at 45-60 degrees.

FW-133



## Mesh Arm Sling

Made from breathable mesh fabric. Has wide, comfortable and adjustable straps. Is compatible with both right and left sides.

FW-134



## Clavicle Bandage

Provides stabilization in clavicle and sternum-clavicle fractures. It has compression and pain relieving effects. Its 8-shape design and adjustable velcro straps make it easy to put on and take off.

FW-135



## Arm Sling (Lux)

Made from fabric that keeps warm. Has wide, comfortable and adjustable straps. Thanks to the waist belt support, there is no excessive load on the neck and other shoulder. Is compatible with both right and left sides.



FW-136

## Shoulder Support

Made from neoprene material. Provides warmth and support to the elbow joint. With its pressure effect, provides an antalgic effect with its pressure effect.



FW-137

## Shoulder Support (Standard)

Made from 3mm neoprene material. Applies compression. Prevents edema and reduces pain by supporting the shoulder from the armpit area.



FW-138

## Velpau Bandage

Designed to keep the arm and shoulder fixed in the desired position in cases of humerus fractures and shoulder dislocation. Has adjustable straps and is fixed with a chest strap. It is practical and designed to be easy to use.



FW-139

## Elbow Support

Made from neoprene material. Provides warmth and support to the elbow joint. Provides an antalgic effect with its pressure effect.

FW-140



## Elbow Support (Standard)

Made from neoprene material. Provides warmth and support to the elbow joint. Provides an antalgic effect with its pressure effect.

FW-141



## Epicondylitis Bandage

It bypasses the elbow and epicondyle by applying pressure to the tendons and reducing the load on the epicondylitis. It can be adjusted according to the person thanks to the band on it.

FW-142



## Epicondylitis Bandage With Silicon Pad

Has a silicon pad that provides local compression and a gel pad that allows for cold therapy use. Bypasses the elbow and epicondyle by applying pressure to the tendons with its designed three-point pad and reduces the load on epicondylitis. Can be adjusted according to the person with the band on it.

FW-143



## Elbow Support With Pad

Made of 3 mm neoprene. With its soft pad, it provides high protection in cases such as olecranon bursitis and sulcus ulnar syndrome. Provides support and protection for athletes.

FW-144



## Epicondylitis Bandage With Gel Pad

Has a gel pad that provides local compression and a gel pad that allows for cold therapy use. Bypasses the elbow and epicondyle by applying pressure to the tendons with its designed three-point pad and reduces the load on epicondylitis. Gel pad provides cold treatment. Can be adjusted according to the person with the band on it.

FW-145



## Tennis Elbow Support

Made from neoprene fabric. Provides support to the elbow joint. With the adjustable band at the bottom, it creates a compression effect on the muscle group connected to the epicondyle and bypasses the epicondyle in the elbow joint.

FW-146



## Elbow Immobilizer

Made from neoprene fabric. Prevents sweating thanks to the laminated towel on its inner surface. Keeps the elbow joint in extension with 3 steel supports.

FW-147



## Cubital Tunnel Elbow Support

Used in cubital tunnel syndrome and ulnar nerve damage. The 40-degree inclined steel support prevents bending of the elbow and stretching of the nerve, so that the nerve rests in its maximum comfortable position without being trapped.

FW-148



## Angle Adjusted Brace

Designed in accordance with the legs anatomical shape. There are aluminum supports on the medial and lateral sides of the orthosis. There is a hinged joint and ROM degree adjustable device in the knee joint section. Flexion/extension angle is adjusted with the angle adjustment mechanism. The aluminum bar and the surfaces in contact with the skin are covered with terry cloth.



FW-149

## Camel Hair Wool Knee Pad

Supports joint function. Due to its superior knitting technology, it wraps the knee tightly and locally from four sides, provides heat increase and creates a proprioceptive effect. Its knitted spacing and elastic material provide effective support to weak or tight knee joints. Provides protection against new damages that may occur during the continuation of active mobility. Used for arthritis and rheumatism pain.



FW-150

## Simple Knee Support

Made from neoprene material. Provides analgesic effect with compression. Supports the weak muscles around the knee. Provides stabilization in the knee joint.



FW-151



## Immobilizer (45-55-65cm)

It is made from sponge fabric for patient comfort. There are removable medial, lateral and posterior aluminum supports that provide immobilization of the knee joint in extension. Produced in 3 sizes according to leg length.



FW-152

## Simple Knee Support (Standard)

Made from neoprene material. Provides analgesic effect with compression. Supports the weak muscles around the knee. Provides stabilization in the knee joint. It is comfortable to wear with its Velcro system and can be adjusted to any size.



FW-153

## Short Patella Supported Knee Brace (Standard)

Made from neoprene material. The area around the patella is open and supported by a pad. Provides analgesic effect with compression. Supports the weak muscles around the knee. Provides stabilization in the knee joint. It is comfortable to wear with its Velcro system and can be adjusted to any size.



FW-154

## Short Knee Brace With Patella And Ligament Support (Standard)

Made from neoprene material. The area around the patella is open and supported by a pad. In medial and lateral, it is supported by four spiral springs. Provides stabilization in the knee joint with flexible underwires on the side. Supports the ligaments and relieves the load on the knee joint. Provides analgesic effect with compression. Supports the weak muscles around the knee.

FW-155



## Hinged Knee Brace

Made from neoprene material. The area around the patella is open. It is suitable for use after serious surgeries and injuries of the knee. Provides knee joint stabilization with its removable hinged aluminum support on both sides. Provides analgesic effect with compression. Supports the weak muscles around the knee.

FW-156



## Hinged Knee Brace (Standard)

Manufactured of neoprene material. The area around the patella is open. It is suitable for use after serious surgeries and injuries of the knee. With removable hinged aluminum support on both sides, it allows controlled motion and prevent knee hyperextension. Provides analgesic effect with compression.

FW-157



## Hinged Knee Support (Open Front)

Manufactured of neoprene material. The area around the patella is open. Since it can be fully opened from the front, it provides comfortable use for the elderly and overweight people who have recently undergone surgery. Velcro fasteners on the upper and lower calves increase stabilization. It is suitable for use after serious surgeries and injuries of the knee.

FW-158



## Hinged Knee Support (Open Top)

Manufactured of neoprene material. The area around the patella is open. It is suitable for use after serious surgeries and injuries of the knee. With removable hinged aluminum support on both sides, it allows controlled motion and prevent knee hyperextension. Provides analgesic effect with compression.

FW-159



## ACL Supported Knee Brace

Made from neoprene material. The area around the patella is open. Flexible underwires and cross bands on both sides provide stabilization of the cruciate ligaments, joints, and patella.

FW-160



## Patella Supported Knee Brace

Made from neoprene material. The area around the patella is open and supported by a pad. Provides analgesic effect with compression. Supports the weak muscles around the knee. Provides stabilization in the knee joint.

FW-161



## Patella Supported Knee Brace (Standard)

Made from neoprene material. The area around the patella is open and supported by a pad. Provides analgesic effect with compression. Supports the weak muscles around the knee. Provides stabilization in the knee joint. It is comfortable to wear with its Velcro system and can be adjusted to any size.

FW-162



## Knee Brace With Patella And Ligament Support

Made from neoprene material. The area around the patella is open and supported by a pad. In medial and lateral, it is supported by four spiral springs. Provides stabilization in the knee joint with flexible underwires on the side. Supports the ligaments and relieves the load on the knee joint. Provides analgesic effect with compression. Supports the weak muscles around the knee. Provides stabilization in the knee joint.

FW-163



## Patella And Ligament Supported Open Top (Flexible Underwire)

Made from neoprene material. The area around the patella is open and supported by a pad. It is open at the top which provides comfortable use for elderly and overweight patients. In medial and lateral, it is supported by four spiral springs. Provides stabilization in the knee joint with flexible underwires on the side. Supports the ligaments and relieves the load on the knee joint. Provides analgesic effect with compression. Provides stabilization in the knee joint and supports the weak muscles around the knee.

FW-164



## Patellar Tendon Strap (Standard)

Designed and manufactured from neoprene material. Has a pad that supports the patella tendon. Provides analgesic effect with compression. It is suitable for Osgood-Schlatter, patellar chondromalacia and patellar tendon disease.


FW-165



## Patellar Tendon Knee Brace

Made from neoprene material. Has a pad that supports the patella tendon. Provides analgesic effect with compression. It is suitable for Osgood-Schlatter, patellar chondromalacia and patellar tendon disease.

FW-166



## Knee Braces For Athletes

It is made of neoprene material. Patella and its surroundings are made of EVA material and supported with pads. Provides constant compression and warmth while protecting the knee against impacts. Maintains body temperature, accelerates blood circulation and thus accelerates the healing process. Protects the patella and soft tissues during sports activities.

FW-167



## Long Knee Brace With Patella And Ligament Support

Made from neoprene material. The area around the patella is open and supported by a pad. In medial and lateral, it is supported by four spiral springs. Provides stabilization in the knee joint with flexible underwires on the side. Supports the ligaments and relieves the load on the knee joint. Provides analgesic effect with compression. Provides stabilization in the knee joint and supports the weak muscles around the knee.

FW-168



## Long Knee Brace With Patella And Ligament Support (Standard)

Made from neoprene material. The area around the patella is open and supported by a pad. In medial and lateral, it is supported by four spiral springs. Provides stabilization in the knee joint with flexible underwires on the side. Supports the ligaments and relieves the load on the knee joint. Provides analgesic effect with compression. Provides stabilization in the knee joint and supports the weak muscles around the knee.

FW-169



## Camel Hair Wool Knee Pad

Supports joint function. Due to its superior knitting technology, it wraps the knee tightly and locally from four sides, provides heat increase and creates a proprioceptive effect. Its knitted spacing and elastic material provide effective support to weak or tight knee joints. Provides protection against new damages that may occur during the continuation of active mobility. Used for arthritis and rheumatism pain.

FW-170



## Lace-Up Wrist Rest Splint With Thumb Support (Standard)

Provides thumb immobilization and keeps the wrist in a neutral position. Used to restrict thumb and wrist movement in de Quervain and carpal tunnel syndrome (stabilization of the thumb and wrist joint). It is easy to put on and provide stabilization with its designed laces.

FW-171



## Lace-Up Wrist Splint

Provides immobilization of the wrist in a neutral position. Used in wrist immobilization and carpal tunnel syndrome. Due to its designed laces, it is easy to put on and provide stabilization.

FW-172



## De Quervain Splint

There is an aluminum that keeps the thumb in a neutral position. The aluminum part is extended above the wrist. Used to treat De Quervain syndrome and restrict movement of the thumb and wrist. Provides stabilization of the thumb and wrist joint.

FW-173



## De Quervain Splint (Standard)

There is an aluminum that keeps the thumb in a neutral position. The aluminum part is extended above the wrist. Used to treat De Quervain syndrome and restrict movement of the thumb and wrist. Provides stabilization of the thumb and wrist joint.

FW-174



## Wrist Splint With Thumb Support

Ensures thumb immobilization and the wrist in a neutral position. Used to restrict thumb and wrist movement in de Quervain and carpal tunnel syndrome (stabilization of the thumb and wrist joint).



FW-175

## Wrist Rest Splint With Thumb Support (Double Sided Mesh)

Used for thumb and wrist immobilization after plaster casts and surgery in cases of De Quervain and carpal tunnel disorders. Has bandage around the wrist that increases stabilization. Has two underwire channels in right and left. Metal underwire is placed in one of these channels according to the hand on which the product will be applied. Underwire always falls on the inner part of the hand on which the product is applied.



FW-176

## Wrist Rest Splint With Thumb Support (Standard)

Ensures thumb immobilization and the wrist in a neutral position. Used to restrict thumb and wrist movement in de Quervain and carpal tunnel syndrome (stabilization of the thumb and wrist joint).



FW-177

## Simple Wristband

Made from neoprene material. Provides better protection by wrapping the wrist in the form of a bandage. Helps rehabilitation by preserving heat in wrist injuries, sprains and traumas. Prevents injuries and dislocations that may occur during sports and physical activities.



FW-178

## Boxer Fracture Splint (Standard)

Used for boxer fractures and contusions. Provides stabilization of the 4th and 5th fingers and wrist.



FW-179

## Wrist Rest Splint

Provides immobilization of the wrist in a neutral position. Used in wrist immobilization and carpal tunnel syndrome.



FW-180

## Wrist Rest Splint(Standard)

Provides immobilization of the wrist in a neutral position. Used in wrist immobilization and carpal tunnel syndrome.



FW-181

## Bidirectional Hand Wrist Rest Splint (Standard)

Provides immobilization of the wrist in a neutral position. Used for wrist immobilization and carpal tunnel syndrome. Can be used for both hands.



FW-182



## Thermoplastic Anti-Spasticity Splint

Made from thermoplastic material. Has spaced finger channels. In spastic cases, it keeps the hand in a neutral position by preventing flexion of the hand and fingers.



FW-183

## Thermoplastic Dynamic Hand-Finger Splint

Made from thermoplastic material. Passively stabilizes the wrist and fingers in extension. Allows active finger flexion. Allows controlled movement in the early period to prevent the formation of adhesions after tendon repair in extensor tendon injuries.



FW-184

## Thermoplastic Hand Immobilisation Splint

Made from thermoplastic material and is lightweight. The interior is lined with soft quality fabric. Keeps the wrist at rest in a neutral position.



FW-185

## Thermoplastic Kleinert Splint

Made from thermoplastic material. Passively stabilizes the wrist and fingers in flexion. Active finger allows extension. Allows controlled movement in the early period in order to prevent the formation of adhesions after tendon repair in flexor tendon injuries.



FW-186

## Nelson Collar With Chin Support

Adjusted to anatomical neck shape. It has holes that help skin to breathe. It is made of 12mm anti-allergic plastazote material. Has special support for the jaw. Prevents flexion and keeps the neck in extension. Provides semi-rigid support.

FW-187



## Minerva Collar

A corset that extends to the lower costal margin. Provides good control of all cervical motion. Includes adjustable height. The dorsal part allows the patient to sleep on his/her back without being disturbed. Provides more movement restriction with respect to Philadelphia.

FW-188



## Nelson Collar

Adjusted to anatomical neck shape. Made of 12mm anti-allergic plastazote material. Has holes that help skin to breathe. Restricts neck movements and supports muscles. Provides semi-rigid support.

FW-189



## Philadelphia Collar

Provides full stability during emergency transport of patients with neck fractures. Consists of occiput and mandibular parts. Properly limits neck movements.

FW-190



## Soft Cervical Collar

Adjusted to anatomical neck shape. Produced from 45 density hard sponge material. The fabric part is made of Bamboo, which is soft and has the feature of keeping the skin dry. Provided in two different heights (8 cm and 10 cm). Can be washed due to having the extra cover on it. Has an adjustable Velcro fastener on the back. Keeps the neck warm, soft neck collar that restricts neck movements and supports muscles.

FW-191



## Philadelphia Tracheostomy Collar

Consists of occiput and mandibular parts. Properly limits neck movements. Suitable for tracheostomy.

FW-192



## Adjustable Height Cervical Collar With Chin Rest

Designed according to Neck Ergonomics. Due to its breathable materials, it keeps your neck dry and odour-free. Due to its layer by layer and special Velcro design it offers adjustable height and circumference for neck. Has special support for the jaw. Prevents flexion and keeps the neck in extension.

FW-193



## Height Adjustable Cervical Neck Collar

Designed according to Neck Ergonomics. Due to its breathable materials, it keeps your neck dry and odour-free. Due to its layer by layer and special Velcro design it offers adjustable height and circumference for neck. Restricts neck movements and supports neck muscles. Provides semi-rigid support.

FW-194



## Plaster Shoes

Prevents the plaster base from getting dirty and slipping. Protects the plaster against mechanical damage such as wear and tear. It is a light slipper. It can be easily worn on the cast with elastic bands attached to the velcro straps.

FW-195



## Figure-8 Ankle Support

Made from neoprene material. It increases stabilization with its eight-figure wrapping tapes. Provides analgesic effect with compression. Supports weak muscles and ligaments around the foot and ankle. Provides stabilization in the ankle joint.

FW-196



## Lace Up Ankle Brace

Designed Underwires on the medial and lateral parts of the ankle increase its stability. The designed lace-up system provides better coverage of the foot and ankle and increases stability. Provides analgesic effect with compression. It provides stabilization in the ankle joint by supporting the weak muscles and ligaments around the foot and ankle.

FW-197



## Cross Band Ankle Brace

Made from neoprene material. Increases stabilization with its eight-figure wrapping tapes. Provides analgesic effect with compression. Supports weak muscles and ligaments around the foot and ankle. Provides stabilization in the ankle joint.

FW-198



## Cross Band Ankle Brace (Standard)

Made from neoprene material. Increases stabilization with its eight-figure wrapping tapes. Provides analgesic effect with compression. Supports weak muscles and ligaments around the foot and ankle. Provides stabilization in the ankle joint.

FW-199



## Dorsiflexion Band (Foot Up)

Made from neoprene material. Used in foot drop problems. Keeps the ankle in dorsiflexion at the desired tension.

FW-200



## Elastic Hallux Valgus Day Splint

It is a practical and comfortable splinter for daily use. Can be used inside wide shoes or slippers. Used for conservative treatment of hallux valgus and after surgery. Prevents the progressive deformity by pushing the thumb out.

FW-201



## Foot Guard AFO

Made from thermoplastic material. Covered with removable and washable terry cloth. The T apparatus ensures that the ankle is fixed at 90 degrees and prevents internal/external rotations.

FW-202



## Hallux Valgus Shoes

Provides early mobilization after forefoot surgeries such as hallux valgus. It is also used in diabetic patients to reduce the load on the forefoot in non-healing wounds on the forefoot. Its sole is made of rubber material to prevent slipping. Its upper surface is leather. The slope of the sole reduces the load on the forefoot to a negligible level. Reduces foot edema and shortens the treatment process. The velcro on the front part of product is designed for easy adjustment.

FW-203



## Hallux Valgus Night Splint

Used for conservative treatment of hallux valgus and after surgery. Suitable for night use. It is a therapeutic support that helps heal bunions during sleep. Protects and supports the soft tissue and joint structure in the area by pushing the thumb laterally. Prevents the progression of deformities in the joint and the formation of bunion. Helps the healing of the bunion lesion.

FW-204



## Ankle Brace With Air Cushioning

A product designed with outer plastic part and inside air cushions. Since the inside of the airbag is closed, it will wrap the ankle well according to its anatomy. Provides ankle stability. Due to its all-round airbag, the applied compression can be adjusted as desired. The air pad can be easily inflated using a special mini pump.

FW-205



## Gel Pad Ankle Brace

The outer surface is made of plastic material. It is a foot orthosis the inner surface of which consists of gel pads for hot or cold treatment. The gel pad can be frozen for cold treatment. The gel pad also has the feature of being heated in hot water, microwave or oven for hot treatment.

FW-206



## Achilles Boot (Short)

Offers both a functional and optimal solution for the post-surgical treatment of Achilles ruptures. The Achilles boot provides immobility, protection and comfort. The base is made of hard plastic, the outer part is made of plastic material and the inside is lined with a terry pad. Due to its special sole system, it provides physiological walking. Supports easy walking due to the cradle system design of the sole.



FW-207

## Ankle Brace With Ligament Support

Made from neoprene material. Additionally designed Plastic ligaments for malleol support, increase the stability of the ankle. Provides analgesic effect with compression. Supports weak muscles and ligaments around the foot and ankle. Provides stabilization in the ankle joint.



FW-208

## Ankle Brace With Ligament Support

Made from neoprene material. Additionally designed Plastic ligaments for malleol support, increase the stability of the ankle. Provides analgesic effect with compression. Supports weak muscles and ligaments around the foot and ankle. Provides stabilization in the ankle joint.



FW-209

## Malleol Supported Ankle Brace

Made from neoprene material. Increases stabilization with its elastic bands wrapped in 8-figure. Has pads on both sides that support the inner and outer malleoli. Provides analgesic effect with compression. Provides stabilization in the ankle joint by supporting the weak muscles and ligaments around the foot and ankle.



FW-210

## Malleol Supported Ankle Brace (Standard)

Made from neoprene material. Increases stabilization with its elastic bands wrapped in 8-figure. Has pads on both sides that support the inner and outer malleoli. Provides analgesic effect with compression. Provides stabilization in the ankle joint by supporting the weak muscles and ligaments around the foot and ankle.

FW-211



## Padded Ankle Brace

Designed according to the joint shape, providing internal and external stabilization of the ankle without hindering natural joint movement. Its outer surface is made of plastic material. Because of its sponge pads on the inner surface, it can be used easily in mild and moderate sprains.

FW-212



## Reflex AFO (Right-Left)

It is an ankle orthosis made of plastic material. Provides static dorsiflexion assistance and lateral stability for the ankle/foot area. Velcro fastening that covers the foot and ankle provides comfortable application. For use in shoes, it can be cut according to the inner shape of the shoe and used easily.

FW-213



## Rom Walker (Short)

Provides protected freedom of movement after trauma or surgery. The base is made of hard plastic material and the inside is lined with a terry pad. Due to the angle adjustable mechanism, desired flexion and extension adjustments can be made. The range of motion is adjustable in 0-40 degrees of planta flexion and 0-30 degrees of dorsiflexion. Helps early mobilization and provides natural and comfortable walking with its special sole design. The sole of the foot is wide to accommodate edema or bandages. Inner/outer cushioned sole absorbs shock to the foot and ankle during heel strike.

FW-214



## Rom Walker (Long)

Provides protected freedom of movement after trauma or surgery. The base is made of hard plastic material and the inside is lined with a terry pad. Due to the angle adjustable mechanism, desired flexion and extension adjustments can be made. The range of motion is adjustable in 0-40 degrees of planta flexion and 0-30 degrees of dorsiflexion. Helps early mobilization and provides natural and comfortable walking with its special sole design. The sole of the foot is wide to accommodate edema or bandages. Inner/outer cushioned sole absorbs shock to the foot and ankle during heel strike.

FW-215



## Achilles Boot (Long)

Offers both a functional and optimal solution for the post-surgical treatment of Achilles ruptures. The Achilles boot provides immobility, protection and comfort. The base is made of hard plastic, the outer part is made of plastic material and the inside is lined with a terry pad. Due to its special sole system, it provides physiological walking. Supports easy walking due to the cradle system design of the sole.

FW-216





## Silicone Callus Pads (Standard)

Prevents the formation of calluses caused by friction or pressure. Relieves the pain of existing calluses. Accelerates the healing process by acting as a buffer on the calloused area. Due to its strong adhesive, it remains fixed in place for a long time.



FW-217

## Silicone Insoles With Anatomical Support

Made of special silicone material that does not cause allergies and does not harm the skin. Used as a pain absorber in cases of Fat Pad, Hyperkeratosis and Planter Fasciitis. Provides anatomical support and allows evenly distribution of the pressure on the sole. Provides relief from pain caused by disorders such as Pes Cavus. Relieves partial pressure and tension in the feet caused by standing for a long time and posture disorders. The low-density structure of the blue area increases adaptation to the sole of the foot and protects sensitive and bony areas.



FW-218

## Arch Support Cushions

Made of special silicone material that does not cause allergies and does not harm the skin. Reduces sole stress by supporting the arch area in cases of flat feet. Due to its adhesive feature, it remains stable in the shoe.



FW-219

## Arch Support Socks

Eliminate flat foot discomfort in your feet. Allows you to do your daily activities such as walking and running in the most comfortable way. Wearable arch support socks fill your foot space and perform a supportive function.



FW-220

## Arc Supported Fabric Covered Adhesive Silicone

Made of special silicone material that does not cause allergies and does not harm the skin. Reduces sole stress by supporting the arch area in cases of flat feet. Due to its adhesive feature, it remains stable in the shoe.



FW-221

## Arc Supported Silicone Insoles

Made of special silicone material that does not cause allergies and does not harm the skin. Ensures balanced distribution of body load on the foot. Reduces the load and stress on the joints by absorbing reaction forces. Unlike other insoles, it has silicone support in the arch area.



FW-222

## Splint Bunion Protector

Made with gel pads and elastic fabric. It relieves discomfort caused by metatarsalgia on the sole of the foot and provides relief by massaging the bunion area. Its design ensures that the gel protector always remains in the correct position. It adapts to all foot sizes with its elastic and soft fabric. It can be used easily in all types of shoes.



FW-223

## Shoe Heel Striking Prevention Gel Pad

It is compatible with every shoe model. The gel structure adheres to the shoe on its own. It is reusable. It is especially suitable for use in the heel area. It adapts to your shoes with its flexible and thin structure.



FW-224

## Ingrown Nail Removal Apparatus

Helps you solve your painful ingrown nail problem. Can remove your problematic nail from the skin by pulling it from two directions and reduce your pain in a short time.



FW-225

## Five Toes Bunion Protector

It spreads your fingers apart. The foot cushion prevents slipping between the toes and rash. It is dominated and used by dislocations in the finger bones. The toe protector is suitable for users who stand and work all day. Shoes and socks can be worn inside.



FW-226

## Metatarsal Sleeve With Gel Cushion

It will prevent bone growth (hallux valgus deformity) in your foot. A protector that will help heal deformities in your fingers. It will help you get rid of the problems you experience with regular use. You can wear it under shoes and socks.



FW-227

## Invisible Height Increase Insoles

The product is made of quality silicone material. Increase the height by up to 4 cm. Reduces foot heel pain. Good for pain caused by cracks.



FW-228

## Gel Bunion Corrector Toe Protector

Prevents pressure interruptions caused by shoes. It is used in cases of hallux valgus.



FW-229

## Steel Insoles

Supported by a stainless steel bridge and the upper part is made of leather. It is an insole with metatarsal arch support. It is used in the treatment of Pes Planus, Pes Cavus, Metatarsalgia and Epin Kalkanei deformities.



FW-230

## Bunion Spreader In Gel With Double Ring 1 Pc

Made of high-quality silicone material. Separates the first and second toes with its roller part, prevents the progression of hallux valgus deformity and relieves pain in the area with its bunion protective part. Due to its soft structure, it does not cause any discomfort. Suitable for both right and left feet.



FW-231

## Perforated Heel Protector

Made from qualified silicone material for dynamic energy support and cushioning. With its excellent shock absorption capacity, it relieves foot pain and protects your feet from injury. Protects the soles of the feet from cracking and moisture. Due to its self-adhesive structure, it remains fixed inside the shoe.



FW-232

## Silicone Exercise Ball

Made of 100% silicone materials. Strengthens wrists and fingers for sports purposes. Their hardness varies depending on the colors. You can exercise by squeezing the ball in your palm and releasing it. Hardness levels: Red: Hard, Blue: Medium, Green: Soft



FW-233

## Epin Heel Cup (Polyurethane)

Made of polyurethane material. Absorbs and reduces impacts on the waist, knee and ankle joints. Used in cases of heel spurs.



FW-234

## Epin Heel Cover – Fabric Covered Adhesive Silicone


Made of special silicone material that does not cause allergies and does not harm the skin. Can be used in painful conditions of the heel and in cases of heel spurs. Due to its adhesive feature, it remains stable in the shoe.



FW-235

## Hallux Valgus Socks

Made with gel pads and elastic fabric. It relieves discomfort caused by metatarsalgia on the sole of the foot and provides relief by massaging the bunion area. Its design ensures that the gel protector always remains in the correct position. It adapts to all foot sizes with its elastic and soft fabric. It can be used easily in all types of shoes.



FW-236

## Height Increase Shoe Insoles With Air Cushion

It is easy to use, soft and comfortable. Reduces the pressure on the knee joints while walking or running. Fits easily with all types of shoes. Because of its special structure its height can be adjusted. You can use it for many years.



FW-237

## Inner – Outer Wedge Insoles

Made of special silicone material that does not cause allergies and does not harm the skin. Used to provide balance on the inside and outside feet.



FW-238

## Ingrowing Toenail Protector

Narrow, pointed-toe shoes cause ingrown toenails. By using the Ingrown Nail Protector, you can put an end to your bad days due to ingrown nails. You can put the gel cap on your finger to prevent the shoe from applying pressure. Thanks to its special structure, it takes the shape of your finger and provides comfort all day long.



FW-239

## Hallux Valgus Splints

Made from high quality, durable medical grade silicone material. Can be used while doing yoga or even jogging. Can be used in wide-toed athletic shoes. The toe straightener is perfect for bunion, hallux valgus, Hammer toes, claw toes and other foot problems.



FW-240

## Anti-Slip Front Insoles

Made from high quality, durable medical grade silicone material. Has adhesive on the bottom and does not slip in shoes. Can be washed and reused after use.



FW-241

## Anti Allergic Soft Silicone Calf Pads

Made of soft, permanent, odorless, featureless, skin-friendly silicone. It visually hides your slimming and provides the perfect look on your legs.



FW-242

## Fabric Covered Flip Flop Reel

It is comfortable to use. Thanks to its wearable fabric structure, it provides a more convenient and comfortable use while standing compared to other silicone finger rollers.



FW-243

## Fabric Covered Silicone Tube Finger Protector

Helps with overall toe protection. With its flexible knitted fabric-covered structure, it helps relieve pain caused by calluses, nail problems, dry cuticles and abrasions. It can be easily cut to fit any desired size.



FW-244

## Fabric Covered Silicone Insoles

It is made of special silicone material that does not cause allergies and does not harm the skin. In addition to its silicone insoles, the fact that it is covered with fabric makes its use convenient and comfortable.



FW-245

## Cork Insoles

The lower layer of cork is covered with special leather. Provides longitudinal and metatarsal support.



FW-246

## Cork Insoles (With Epine)

It is a cork orthopedic insole designed according to the anatomical form of the foot. The surface in contact with the foot is covered with high quality leather. Epine support reduces the load on the heel. Its design, which facilitates natural walking, allows the body load to be distributed evenly on the sole of the foot. It reduces the load on the joints. It can be easily used in shoes with its non-slip structure on the bottom surface.



FW-247

## Massage Insoles

Compressed gel parts provide a massage effect and provide comfortable walking. The thin gel sole suitable for all types of shoes can be used easily without making your shoes smaller. By reducing the pressure on the foot, it relieves the knee, ankle and back pain.



FW-248

## Silicone 5-Toe Spacer

Since it has a soft structure, it will not harm your toes at all. It is a high-quality product produced for deformities in the fingers. It can also be used while applying nail polish. It is produced to correct deformities in the fingers. There are two pieces, Right and Left foot.



FW-249

## Silicone Toe Separator

Made of special silicone material that does not cause allergies and does not harm the skin. Prevents the progression of foot deformation and helps eliminate deformations. Allows overlapping fingers to return to their normal position.



FW-250

## Gel Metatarsal Pad With Toe Spreader

It balances the load distribution on the metatarsal heads. Resolves muscle spasm and eliminates the pain caused by muscle contraction. Forms the metatarsal arch of the foot.



FW-251

## Gel Metatarsal Pad With Toe Ring

It balances the load distribution on the metatarsal heads. Resolves muscle spasm and eliminates the pain caused by muscle contraction. Forms the metatarsal arch of the foot.



FW-252

## Anti-Decubitus Silicone Gel Cushion

Prevents and treats bedbugs or pressure wasps. It can also treat burn patients, multiple sclerosis patients, Lou Gehrig's diseases, etc. can treat. They are used to provide comfort and relaxation.



FW-253

## Silicone Arch Support Sock

Eliminate flat foot discomfort in your feet. Allows you to do your daily activities such as walking and running in the most comfortable way. Wearable arch support socks fill your foot space and perform a supportive function.



FW-254

## Silicone Hammer Finger Cushion

It protects the toes and prevents them from curling. It helps in the treatment of hammer finger.



FW-255

## Silicone Epin Heel Cover

Made of special silicone material that does not cause allergies and does not harm the skin. Reduces Achilles tendon tension by acting as a cushion on the heel. Used to relieve the discomfort caused by the bone protrusion in the heel (Epin Calcaneie/heel spur). The low-intensity blue area is designed to harbor heel spurs.



FW-256

## Silicone Gel Insoles

Made of durable, lightweight, soft and shock absorbing material. Is non-slip and allows the feet to breathe. It is quick-drying and sweat-absorbing.



FW-257

## Silicone Callus Toe Ring

Can be used for all kinds of problems related to toes. Used as a preventative in diabetic patients.



FW-258

## Silicone Callus Toe Sleeve

It is used for calluses, ingrown toenails and similar problems related to toes. Diabetes plays a protective role in patients.



FW-259

## Silicone Booties Socks

Used for heel pain, heel spurs, heel cracks and callus problems. It is worn on the feet like socks. Cream or vaseline can be applied inside before wearing. If cream or vaseline is not applied, it should be powdered. Can be used while wearing shoes.



FW-260

## Silicone Wound Dressing Sheet

Produced in 3 different sizes. (10X10cm-15X15cm-15X30cm)



FW-261

## Silicone Insoles


Made of special silicone material that does not cause allergies and does not harm the skin. Provides anatomical support and allows the pressure on the sole to be distributed evenly. Has a relaxing and comfortable effect. Used for pes planus, fingertip deformation and sole pains. Can be used easily while walking, doing sports and working.



FW-262

## Gel Metatarsal Pad

It is produced with special silicone that does not cause any damage or allergies. It is durable and suitable for long-term use. Metatarsal helps relieve pain by absorbing particle pressure. Due to the adhesive Szeliga, it remains stable in the shoe.



FW-263

## Silicone Half Insoles

Made of special silicone material that does not cause allergies and does not harm the skin. Relief the pains caused by disorders such as Pes Cavus. Relieves partial pressure and tension in the feet caused by standing for a long time and posture disorders. Due to its special design, it protects sensitive and bony areas by increasing adaptation to the sole of the foot.



FW-264

## Starflex Insole

Produced from special thermoplastic. It is a leather-covered insole with longitudinal and metatarsal support on top and bottom. Reduces the load and stress on the joints by ensuring the regular development of the standing body graph.



FW-265

## Starflex Insole (With Epine)

Specially made of thermoplastic material. It is an insole with silicone and metatarsal support in the heel part. The use of insoles eases the load on the joints and reduces pain.



FW-266

## Silicone Exercise Ball

and relaxation exercises. In addition to finger exercises, it can be used for relaxation and stress control. Suitable for use by children and adults.



FW-267

## T-Shaped Shoe Heel Striking Prevention Gel Pad

The silicone on the inside of the shoe will not harm your skin and will protect your foot health. It is very practical to use. Simply remove the protective film from the adhesive surface and place it on the inside of your shoe.



FW-268

## Shoe Heel Striking Prevention Pad

The protector prevents your shoe from hitting your foot from behind. Prevents problems such as wounds and blisters on the top of your heels. At the same time, it makes your shoes easier to use by filling the space in the shoe.



FW-269

## 3-Loop Hammer Toe Cushion

It is a protector that will help you get rid of hammer toe discomfort. This protector, which you will wear on your 3 fingers, will keep your fingers in a straight position and support them with its silicone pad. This way, your fingers will not hurt and you will feel relaxed. You can easily use it under clothing such as socks and shoes. The material is flexible and silicone.



FW-270

## UCBL Insoles

It corrects foot deformities and stops the progression of the deformity. It is a type of orthosis that covers the foot and distributes the entire load of the body equally to the foot. It completely grasps the heel and supports the arch of the foot. It is used for flat soles, excessive curved sole, metatarsal collapse, Hallux rigidus (stiff thumb) deformities.



FW-271

## Ring Inger Spreader

Helps heal hammertoe disease. Prevents irritation caused by overlapping fingers. The roller you place between two fingers will not disturb you due to its silicone structure.



FW-272

## Little Finger Roller Bunion Protector

Made of special silicone material that does not cause allergies and does not harm the skin. Due to the bunion protector, you can prevent deformities in your little finger. Takes the shape of your toe and provides comfort all day long. Silicone finger protector can be washed and reused. Due to the special silicone structure, no size or number is required, it fits all foot sizes.

FW-273



## Splint Roller Bunion Protector

Due to its high-quality silicone structure, it will not disturb your feet. In addition to helping your toes stand correctly with its roller structure, it will help your feet heal quickly by preventing deformation. Socks can be worn under shoes.

FW-274



## Adhesive Gel Heel Pad (Men)

Provides evenly distribution of the pressure to the foot. Minimize the pain and impacts in the heels. The gel used in the blue areas of the insole is softer than the other areas. The blue areas are specially made for heel spur discomfort. It is especially recommended for those who have to stand for long periods of time during the day, athletes and those with heel spur disease.

FW-275



## Adhesive Gel Heel Pad (Women)

Provides evenly distribution of the pressure to the foot. Minimize the pain and impacts in the heels. The gel used in the red areas of the insole is softer than the other areas. The red areas are specially made for heel spur discomfort. It is especially recommended for those who have to stand for long periods of time during the day, athletes and those with heel spur disease. Due to its self-adhesive structure, it remains fixed inside the shoe.

FW-276







## Baseball Splint

The inside is covered with plastazote.  
Aluminum support provides immobilization  
of the finger in extension.

FW-277



## Butonier Splint

It is a finger splint made of thermoplastic.  
Keeps the finger in extension at 3 points.

FW-278



## Hammer Finger Splint (Mallet Finger)

Used in cases of fingertip extensor  
tendon injury. Keeps the tip of the  
finger in extension.


FW-279



## Four Winged Finger Splint

Used to fix the finger in the desired position in case  
of finger fractures and injuries. Has four winged  
finger rests. The inner surface is coated with soft  
plastazote and is made of flexible metal material  
that can be easily shaped.

FW-280



## Two-Winged Finger Stabilization Splint

The inside is covered with plastozone. Aluminum support provides immobilization of the finger in extension.



FW-281

## Frog Splint

The inside is covered with plastozone. Aluminum support provides immobilization of the finger in extension. Frog splints keep the distal interphalangeal joint in extension and the proximal interphalangeal joint in an anatomically appropriate flexion position.



FW-282

## Middle Finger Fracture Splint (Standard)

Aluminum support provides immobilization of the finger in extension.



FW-283

## Knee High Varicose Stocking (Open/Ccl I)

Provides a mild superficial effect. Used in mild varicose veins without edema. Applied for mild varicose veins that occur during and after pregnancy. Used as a healing aid after venous surgery. Used in all kinds of travel.  
Class 1 pressure sock, 15-20 mmHg.

## Knee High Varicose Stocking (Open Toe/Ccl II)

Provides moderate superficial effect. Used in edematous and prominent varicose veins. Applied in venous insufficiency caused by varicosities and after healing of simple ulcers. Used for varicose veins that occur during and after pregnancy. Used after sclerotherapy and stripping.  
Class 2 pressure sock, 20-30 mmHg.

## Knee High Varicose Stocking (Open Toe/Ccl III)

Provides superficial and deep effects. Used in edematous and prominent varicose veins. Used in chronic venous insufficiency due to post-traumatic syndrome. Applied in post-traumatic edema and venous ulcerations. Applied for secondary varicosities after stasis dermatitis. Used after sclerotherapy and stripping in white atrophy.  
Class 3 pressure sock, 30-45 mmHg.



## Thigh High Varicose Stocking (Open Toe/Ccl I)

Provides a mild superficial effect. Used in mild varicose veins without edema. Applied for mild varicose veins that occur during and after pregnancy. Used as a healing aid after venous surgery. Used in all kinds of travel.  
Class 1 pressure sock, 15-20 mmHg.

## Thigh High Varicose Stocking (Open Toe/Ccl II)

Provides moderate superficial effect. Used in edematous and prominent varicose veins. Applied in venous insufficiency caused by varicosities and after healing of simple ulcers. Used for varicose veins that occur during and after pregnancy. Used after sclerotherapy and stripping.  
Class 2 pressure sock, 20-30 mmHg.

## Thigh High Varicose Stocking (Open Toe/Ccl III)

Provides superficial and deep effects. Used in edematous and prominent varicose veins. Used in chronic venous insufficiency due to post-traumatic syndrome. Applied in post-traumatic edema and venous ulcerations. Applied for secondary varicosities after stasis dermatitis. Used after sclerotherapy and stripping in white atrophy.  
Class 3 pressure sock, 30-45 mmHg.



FW-285

## Pantyhose Varicose Stocking (Open Toe/Ccl I)

Provides a mild superficial effect. Used in mild varicose veins without edema. Applied for mild varicose veins that occur during and after pregnancy. Used as a healing aid after venous surgery. Used in all kinds of travel.  
Class 1 pressure sock, 15-20 mmHg.

## Pantyhose Varicose Stocking (Open Toe/Ccl II)

Provides moderate superficial effect. Used in edematous and prominent varicose veins. Applied in venous insufficiency caused by varicosities and after healing of simple ulcers. Used for varicose veins that occur during and after pregnancy. Used after sclerotherapy and stripping.  
Class 2 pressure sock, 20-30 mmHg.

## Pantyhose Varicose Stocking (Open Toe/Ccl III)

Provides superficial and deep effects. Used in edematous and prominent varicose veins. Used in chronic venous insufficiency due to post-traumatic syndrome. Applied in post-traumatic edema and venous ulcerations. Applied for secondary varicosities after stasis dermatitis. Used after sclerotherapy and stripping in white atrophy.  
Class 3 pressure sock, 30-45 mmHg.



## Knee High Varicose Stocking (CloseToe/Ccl I)

Provides a mild superficial effect. Used in mild varicose veins without edema. Applied for mild varicose veins that occur during and after pregnancy. Used as a healing aid after venous surgery. Used in all kinds of travel.  
Class 1 pressure sock, 15-20 mmHg.

## Knee High Varicose Stocking (CloseToe/Ccl II)

Provides moderate superficial effect. Used in edematous and prominent varicose veins. Applied in venous insufficiency caused by varicosities and after healing of simple ulcers. Used for varicose veins that occur during and after pregnancy. Used after sclerotherapy and stripping.  
Class 2 pressure sock, 20-30 mmHg.

## Knee High Varicose Stocking (CloseToe/Ccl III)

Provides superficial and deep effects. Used in edematous and prominent varicose veins. Used in chronic venous insufficiency due to post-traumatic syndrome. Applied in post-traumatic edema and venous ulcerations. Applied for secondary varicosities after stasis dermatitis. Used after sclerotherapy and stripping in white atrophy.  
Class 3 pressure sock, 30-45 mmHg.



FW-287

## Thigh High Varicose Stocking (Close Toe/Ccl I)

Provides a mild superficial effect. Used in mild varicose veins without edema. Applied for mild varicose veins that occur during and after pregnancy. Used as a healing aid after venous surgery. Used in all kinds of travel.  
Class 1 pressure sock, 15-20 mmHg.

## Thigh High Varicose Stocking (Close Toe/Ccl II)

Provides moderate superficial effect. Used in edematous and prominent varicose veins. Applied in venous insufficiency caused by varicosities and after healing of simple ulcers. Used for varicose veins that occur during and after pregnancy. Used after sclerotherapy and stripping.  
Class 2 pressure sock, 20-30 mmHg.

## Thigh High Varicose Stocking (Close Toe/Ccl III)

Provides superficial and deep effects. Used in edematous and prominent varicose veins. Used in chronic venous insufficiency due to post-traumatic syndrome. Applied in post-traumatic edema and venous ulcerations. Applied for secondary varicosities after stasis dermatitis. Used after sclerotherapy and stripping in white atrophy.  
Class 3 pressure sock, 30-45 mmHg.





## Pantyhose Varicose Stocking (Close Toe/Ccl I)

Provides a mild superficial effect. Used in mild varicose veins without edema. Applied for mild varicose veins that occur during and after pregnancy. Used as a healing aid after venous surgery. Used in all kinds of travel.  
Class 1 pressure sock, 15-20 mmHg.

## Pantyhose Varicose Stocking (Close Toe/Ccl II)

Provides moderate superficial effect. Used in edematous and prominent varicose veins. Applied in venous insufficiency caused by varicosities and after healing of simple ulcers. Used for varicose veins that occur during and after pregnancy. Used after sclerotherapy and stripping.  
Class 2 pressure sock, 20-30 mmHg.

## Pantyhose Varicose Stocking (Close Toe/Ccl III)

Provides superficial and deep effects. Used in edematous and prominent varicose veins. Used in chronic venous insufficiency due to post-traumatic syndrome. Applied in post-traumatic edema and venous ulcerations. Applied for secondary varicosities after stasis dermatitis. Used after sclerotherapy and stripping in white atrophy.  
Class 3 pressure sock, 30-45 mmHg.



FW-289

## Pregnancy Varicose Socks

Used in cases of varicose vein pain experienced during pregnancy. Provides blood circulation by applying compression to the leg and calf parts during the treatment of high-level varicose veins. It is very useful for varicose veins, which disturb the expectant mother during pregnancy due to the appearance it gives and the pain it causes. Has a relaxing feature and prevents the disease from progressing.



**SOCKS**

**13**

FW-290

## Pediatric Wrist Splint With Thumb Support

Ensures thumb immobilization and the wrist in a neutral position. Used to restrict thumb and wrist movement in de Quervain and carpal tunnel syndrome (stabilization of the thumb and wrist joint).



FW-291

## Pediatric Flat Foot Brace

Designed to be used in mild sprains of the ankle. Provides continuous compression to disperse post-sprain edema. The extra bandage that comes with the product provides extra support to the ankle.



FW-292

## Pediatric Wrist Rest Splint

Provides immobilization of the wrist in a neutral position. Used in wrist immobilization and carpal tunnel syndrome.



FW-293

## Pediatric Mesh Arm Sling

Made from breathable mesh fabric. Has wide, comfortable and adjustable straps. Allows children to fix the arm-shoulder immobilization by adjusting the arm as desired. It is compatible with both right and left points.



FW-294

## Pediatric Immobilizer

Made of flexible fabric for patient comfort. There are medial, lateral and posterior supports that provide immobilization of the knee joints in extension. Manufactured in 4 sizes according to leg length.

FW-295



## Pediatric Clavicle Bandage

Provides stabilization in clavicle and sternum-clavicle fractures. Has compression and pain-relieving effects. Its 8-shape design and adjustable velcro straps make it easy to put on and take off.

FW-296



## Pirate Eye Patch

Pirate Eye Patch is made of high-quality fabric that does not irritate the eyes. Because of its special hollow design, it does not put pressure on the eyes. Can also be used for simple and effective treatments for Amblyopia, Astigmatism, and strabismus. Can be used in both eyes.

FW-297



## Pediatric Ligament Supported Ankle Brace

Made from neoprene material. Plastic ligaments which are additionally designed to support the malleoli increase ankle stability. Provides an analgesic effect with compression. Provides stabilization in the ankle joint by supporting the weak muscles and ligaments around the foot and ankle.

FW-298



## Pediatric Malleol Supported Ankle Brace

Made from neoprene material. Increases stabilization with its elastic bands wrapped in 8-figures. Has pads on both sides that support the inner and outer malleoli. Provides analgesic effect with compression. Provides stabilization in the ankle joint by supporting the weak muscles and ligaments around the foot and ankle.

FW-299



## Pediatric Patella Supported Knee Brace

Made from neoprene material. The area around the patella is open and supported by a pad. Provides analgesic effect with compression. Supports the weak muscles around the knee. Provides stabilization in the knee joint.

FW-300



## Pediatric Knee Brace With Patella And Ligament Support

Made from neoprene material. The area around the patella is open and supported by a pad. Supported by four spiral springs, medial and lateral. Provides joint stabilization with flexible underwires on both sides. Supports the ligaments and relieves the load on the knee joint. Provides an analgesic effect with compression. Supports the weak muscles around the knee. Provides stabilization in the knee joint.

FW-301



## Pediatric Soft Cervical Collar

It is a soft neck collar that restricts neck movements and supports the muscles. Manufactured according to the neck anatomical shape. Produced from 45 density hard sponge material. The fabric part is made of bamboo thread, which is soft and has the ability to keep the skin dry. There is an adjustable velcro at the back. Keeps the neck warm. Has 2 different heights: 8 cm and 10 cm.

FW-302



## Pediatric Velpau Bandage

Designed to keep the arm and shoulder fixed in the desired position in cases of humerus fractures and shoulder dislocation. It is practical, has adjustable straps and is fixed with a chest strap.

FW-303



## Pirate Eye Patch

Pirate Eye Patch is made of high-quality fabric that does not irritate the eyes. Because of its special hollow design, it does not put pressure on the eyes. Can also be used for simple and effective treatments for Amblyopia, Astigmatism, and strabismus. Can be used in both eyes.



FW-304

## Sleep Strip

Made from soft quality fabric that does not irritate the eyes. Is suitable for all head shapes due to its elastic back. Can be used to rest the eyes, protect from light or when suffering from insomnia.



FW-305

## Prominent Ear Bandage

Helps to protect the operated area from external influences after parotid surgery. The doctor and nurse have no problem adjusting the required tightness (pressure). With its special design, there is no slipping in the area where it is applied. Because of the velcro on it, it can be taken off and put on whenever desired. Thus, it can be used easily in long-term treatments.



FW-306

## Anti Snore Chin Strap

Adapts to all facial contours with its velcro strap system. Helps to prevent chin sagging and breath through the nose. Can also be used after aesthetic operations.



FW-307

## Slimming Tight

Made from neoprene material. The zipper on the side provides ease of use during putting on and taking off. In sports activities it increases blood circulation and metabolism by maintaining body temperature. Reduces edema and supports muscles with its compression effect. Increases performance by supporting muscles. Provides an analgesic effect by providing compression. Used for adductor tendinitis and hip muscle injuries.




FW-308





# Aurlina Gel

FW-310



SIZE

60X40X12/8  
55X35X11/7  
60X30X11/7  
50X30X11/7  
45X30X10/6

## ÖZELLİKLER:

ORPEDOX AURLINA JEL Boyun Destekli Visco Yastık yoğunluğu ve duyarlılığı yüksek elastik visco süngerden üretilmiş olup kılıfı için gözenekli pamuklu örme kumaş kullanılmıştır. Özel tasarımı sayesinde boyun bölgesini destekleyerek boşlukta kalmasını engeller ve böylece olası boyun ağrılarını azaltır, rahat bir uyku uyumanıza yardımcı olur. Isı ve basınca göre hareket eden visco elastik (hafızalı sünger) teknolojisi boyun ve omuz kaslarını destekler, sinir hücrelerinin zorlanmasını asgari düzeye indirir.

## FEATURES:

ORPEDOX AURLINA GEL Neck Supported Visco Pillow is made of elastic visco sponge with high density and sensitivity, and a porous cotton knitted fabric is used for its cover. Thanks to its special design, it supports the neck area and prevents it from being in the gap, thus reducing possible neck pain and helping you sleep comfortably. Visco elastic (memory sponge) technology, which acts according to heat and pressure, supports neck and shoulder muscles and minimizes the strain on nerve cells.

## تازي جلا

ORPEDOX AURLINA GEL Neck Supported Visco ينعون صوم تاذ قزوم Visco ةجفنفسا نم متهيو ، ةيساسحو ةيلاع ةفاشك يماسم ينطق جيسن مادختسا افضال غل كويبحم وطف ، صاخلا مديصت لصفب اهعنديو قبقورلا ققطنم معندي ، فسيوجتلا يف دجاوتلا نم قبقورلا مال نم للقي يلاتلابو مونلا ملع كدعاسيو ةلمتحملا حيرم لكشب ةيطاطملا Visco ةينقت يتلا ، (قركاذلا ةجفنفسا)

# Varius

FW-311



60x40x14/12  
52x32x12/10

SIZE

## ÖZELLİKLER:

Çift taraflı tasarımı sayesinde boyun bölgesini destekleyerek boşlukta kalmasını engeller ve böylece olası boyun ağrılarını azaltır, rahat bir uyku uyumanıza yardımcı olur.

Isı ve basınca göre hareket eden visco elastik (hafızalı sünger) teknolojisi boyun ve omuz kaslarını destekler, sinir hücrelerinin zorlanmasını asgari düzeye indirir.

Uyku sırasında yüze uygulanan basınç ile kan dolaşımını düzenleyerek kırıksıklık ve sabah şişkinliklerini azaltır, sağlıklı ve taze bir cilde kavuşmanıza katkı sağlar.

## FEATURES:

Thanks to its double-sided design, it supports the neck area and prevents it from being in the gap, thus reducing possible neck pain and helping you sleep comfortably.

Visco elastic (memory sponge) technology, which acts according to heat and pressure, supports neck and shoulder muscles and minimizes the strain on nerve cells.

By regulating blood circulation with the pressure applied to the face during sleep, it reduces wrinkles and morning puffiness, and contributes to a healthy and fresh skin.

## تأثيره


نیهجولا ىلع همیصرت لصفب بقورلا ققطنم معدی ولف ، غارفلا یف دجاوتلا نم اهنمیو مالاً نم للقی یلاسلابو ، كدعاسیو ؤلمتحملا بقورلا حیدم لكشب جونلا ىلع ؤیطاطملا Visco ؤینقت

یثلا ، (قرئانلا ؤچنفسرا) طغضلاو قرارحلل أقفو لمعت فستكل او بقورلا تالضع معدت ایالخل داہج! نم للقتو ؤیبصللا

ؤیودنلا قروندا میظنت لالخنم انشا هجولا ىلع طغضلاب دىعاجتلا للقی نواف ، جونلا مهاسیو ، یحابصلال خافتنلال او قرشب ىلع لوصولال یف قرضنو ؤیحص

# Varius Jet

FW-312



SIZE

60x40x14/12  
52x32x12/10

## ÖZELLİKLER:

Çift taraflı tasarımı sayesinde boyun bölgesini destekleyerek boşlukta kalmasını engeller ve böylece olası boyun ağrılarını azaltır, rahat bir uyku uyumanıza yardımcı olur.

Isı ve basınca göre hareket eden visco elastik (hafızalı sünger) teknolojisi boyun ve omuz kaslarını destekler, sinir hücrelerinin zorlanmasını asgari düzeye indirir.

Uyku sırasında yüze uygulanan basınç ile kan dolaşımını düzenleyerek kırıksıklık ve sabah şişkinliklerini azaltır, sağlıklı ve taze bir cilde kavuşmanıza katkı sağlar.

## FEATURES:

Thanks to its double-sided design, it supports the neck area and prevents it from being in the gap, thus reducing possible neck pain and helping you sleep comfortably.

Visco elastic (memory sponge) technology, which acts according to heat and pressure, supports neck and shoulder muscles and minimizes the strain on nerve cells.

By regulating blood circulation with the pressure applied to the face during sleep, it reduces wrinkles and morning puffiness, and contributes to a healthy and fresh skin.

## تأثيره

نیه چولای یلع مدمصرت لصفب قبقرلا قطنم مدعی ودف ، غارفلا یف دجاولتال نم ادعنم یو مالاً نم للقی یلاتل ابو ، كدعاس یو قلمت حمل قبقرلا حیدرم لكشرب موزلا یلع قیطاطملا Visco قیقت

یبتلا ، (قركانلا عجنفسرا) ، طغضلاو قرارحلل أقسو لمعت فستكلاو قبقرلا تالضرع مدعت ایالخل داہج! نم للقتو قیصصعلا

قیدمدلا قروندلا دیظنت لال خ نم ءانثا دچولا یلع طغضلاب دیعاجتلا للقی مناف ، موزلا مءاس یو ، حابصلال خافتنالاو قرشب یلع لوصحل یف أقسو لمعت قرضنو قیحص مدعت ، طغضلاو قرارحلل فستكلاو قبقرلا تالضرع ایالخل داہج! نم للقتو قیصصعلا

قیدمدلا قروندلا دیظنت لال خ نم ءانثا دچولا یلع طغضلاب دیعاجتلا للقی مناف ، موزلا مءاس یو ، حابصلال خافتنالاو قرشب یلع لوصحل یف أقسو لمعت قرضنو قیحص مدعت ، طغضلاو قرارحلل فستكلاو قبقرلا تالضرع ایالخل داہج! نم للقتو قیصصعلا

# Aperi

FW-313



60x40x16

SIZE

## ÖZELLİKLER:

ORPEDOX APERİ Yastık nefes alabilen gözenekli yapısı sayesinde hava sirkülasyonunu güçlendirerek baş ve omuz bölgesi ile yastık arasındaki ısı depolanmasını azaltarak kullanıcı için daha ferah bir ortam sunar. Visco elastik (hafızalı sünger) teknolojisi ile üretilmiştir. Isı ve basınca göre hareket eder. Bu sayede boyun ve omuz kaslarını destekleyerek olası ağrı sorunlarını önler. APERİ Yastık elastik yapısı ile başınızın şeklini alarak boynunuza ve omzunuza tam bir destek sağlar. Baş ve sinirlerin zorlanmasını asgari düzeye indirerek, eklem

## FEATURES:

Thanks to its breathable porous structure, ORPEDOX APERİ Pillow provides a more spacious environment for the user by strengthening the air circulation and reducing the heat storage between the head and shoulder area and the pillow. It is produced with visco elastic (memory sponge) technology. It moves according to heat and pressure. In this way, it prevents possible pain problems by supporting the neck and shoulder muscles. APERİ Pillow takes the shape of your head with its elastic structure and provides full support to your neck and shoulder.

## تأثيري

يتماسمها اهلكتيه لضفب رفوت ، سفنثلل لباقلا ORPEDOX APERİ قداسو مدختسملل اعاسرتا رشكأ قئيب ءاولا نارود زيزعت لالخنم نيب قرارحل نيزخت ليلقتو قداسول او فحتكلاو سارلا ققطنم ولسيف قئيقتب هجاتن! متي لمعي .(قرلكاذل جنفسا) قزمل ا هذهبو .طغضلاو قرارحل بسح لكاشم عنمي ناف ، ققيرطلا مدخلالخنم قلمتحملا ملال انيفحتكلاو قبقرلا تالضع كسار لكش ذخات APERİ قداسو مدخل رفوتو نرملا اهلكتيهب نع .كئيفحتكو كفتبقرل لمالكلا يلع طغضلا ليلقت قيرط دعاسي ناف ، باصعالاو سارلا لصافملا مالأ فيفخت يف يرقفلا كدومع يلع ظافحل او ميقتسم طخ يف

# Aperi Jet

FW-314



60x40x16

SIZE

## ÖZELLİKLER:

ORPEDOX APERİ Yastık nefes alabilen gözenekli yapısı sayesinde hava sirkülasyonunu güçlendirerek baş ve omuz bölgesi ile yastık arasındaki ısı depolanmasını azaltarak kullanıcı için daha ferah bir ortam sunar.

Visco elastik (hafızalı sünger) teknolojisi ile üretilmiştir. Isı ve basınca göre hareket eder. Bu sayede boyun ve omuz kaslarını destekleyerek olası ağrı sorunlarını önler. APERİ Yastık elastik yapısı ile başınızın şeklini alarak boynunuza ve omzunuza tam bir destek sağlar. Baş ve sinirlerin zorlanmasını asgari düzeye indirerek, eklem

## FEATURES:

Thanks to its breathable porous structure, ORPEDOX APERİ Pillow provides a more spacious environment for the user by strengthening the air circulation and reducing the heat storage between the head and shoulder area and the pillow.

It is produced with visco elastic (memory sponge) technology. It moves according to heat and pressure. In this way, it prevents possible pain problems by supporting the neck and shoulder muscles. APERİ Pillow takes the shape of your head with its elastic structure and provides full support to your neck and shoulder.

## تأثيري ملأ

يأس ملأ اهل كيه لضفب رفوت ، سفن تل لباقلا مدخ تس ملأ اع اس تا رشكأ قئيب ءاولا نارود زيزعت لال خ نم نيب قراح ل نيزخت ليلقتو قداسول او فحتكلاو سارلا ققطنم ولسيف ءينقتب هجاتنإ متي لمعي . (قوكاذلأ جنفسا) قزملأ مذهبو . طغضلاو قراحلا بسح لكاشم عنمي نإف ، ققيرطلا مدخلال خ نم قلمت حملأ ملأال نيفتكلأو قبقولأ تالضع كسار لكش ذخات APERİ قداسو مدخلأ رفوتو نرملا اهل كيهب نع . كئيفتكو كئتبقرل لمكلا يلع طغضلا ليلقت قيرط دعاسي نإف ، باصعال او سارلا لصافملأ مالأ فيفخت يف يرقفلا كدومع يلع ظافحلأو ميقتسم طخ يف

# Hella



FW-315



60x40x13

**SIZE**

## ÖZELLİKLER:

ORPEDOX HELLA Yastık nefes alabilen gözenekli yapısı sayesinde hava sirkülasyonunu güçlendirerek baş ve omuz bölgesi ile yastık arasındaki ısı depolanmasını azaltarak kullanıcı için daha ferah bir ortam sunar.

Visco elastik (hafızalı sünger) teknolojisi ile üretilmiştir. Isı ve basınca göre hareket eder. Bu sayede boyun ve omuz kaslarını destekleyerek olası ağrı sorunlarını önler. Baş destek bölümünde yer alan kanallar sayesinde hava, ısı ve nem transferlerinin sağlanmasına yardımcı olur.

## FEATURES:

ORPEDOX HELLA Pillow offers a more spacious environment for the user by strengthening air circulation thanks to its breathable porous structure, reducing heat storage between the head and shoulder area and the pillow.


It is produced with visco elastic (memory sponge) technology. It acts according to heat and pressure. In this way, it prevents possible pain problems by supporting the neck and shoulder muscles. Thanks to the channels in the head support section, it helps to provide air, heat and moisture transfers.

## تازي جلا

ORPEDOX HELLA قدايرو رفوت اعاستا رتثكا ةييب زيزعت لال خ نم مدختس ملل اطلشيء لصفب ءاولا نارود امم ، سفن تلل لباقلا يم اس جلا نيب قرار جلا نيزخت نم للقي قداصول او فتكل او سارل ققطنم ولسيف ةينقتب هجاتن ا متي لم عي . (قرلك اذلا جنفسا) قنر جلا هذهبو . طغضل او قرار جلا بسج لكاشم عنمي ناف ، ققير طلا معد لال خ نم قلمت جلا لال ا نيفتكل او قبقرل ا تالضع قنوجوملا تاونقلا لصفب يهف ، سارل ا معد مسق يه ءاولا لقن ريفوت يه دعاست قبوطرل او قرار جلا او

# Hella Jet

FW-316



**SIZE**

60x40x13

## ÖZELLİKLER:

ORPEDOX HELLA jel Yastık nefes alabilen gözenekli yapısı sayesinde hava sirkülasyonunu güçlendirerek baş ve omuz bölgesi ile yastık arasındaki ısı depolanmasını azaltarak kullanıcı için daha ferah bir ortam sunar.

Visco elastik (hafızalı sünger) teknolojisi ile üretilmiştir. Isı ve basınca göre hareket eder. Bu sayede boyun ve omuz kaslarını destekleyerek olası ağrı sorunlarını önler.





## FEATURES:

ORPEDOX HELLA gel Pillow provides a more spacious environment for the user by strengthening the air circulation thanks to its breathable porous structure, reducing the heat storage between the head and shoulder area and the pillow.


It is produced with visco elastic (memory sponge) technology. It moves according to heat and pressure. In this way, it prevents possible pain problems by supporting the neck and shoulder muscles.

## تأثيره

ORPEDOX HELLA جلد قداسو رفوت اعامتار رثكأ ؤئيب ؤيوقت لال خ نم مدختس ملل اهل كيه لصفب ؤاولا نارود امم ، سفن تلل لباقلا يم اس جلا نيب قراحلا نيزخت نم للقي قداصولاو فتكل او سارلا ققطنم ولسيف ؤينقتب هجاتن اءتي (قراكذل جنفسا) قنر جلا قراحلا بسرح كرحتي مناف ، ققيرطلا هذمبو . طغضلا او قلمت جلا ملال لكاشم عنمي قبقرلا تالضع مدع لال خ نم نيفتكل او

COLLECTION	PATTERN	GRM2
<p><b>Bamboo</b> hypoallergenic</p>		<p>300</p>
<p><i>Olive Oil</i></p>		<p>220</p>
<p>COLLAGEN </p>		<p>440</p>


FABRICS



**INCORRECT SLEEPING POSITION**



COLLECTION	PATTERN	GRM2
 Feels so right		300
		240
<b>GEL GROUP</b> FRONT  BACK 		FRONT 180 BACK 240



**CORRECT SLEEPING POSITION**



**Bamboo**  
hypoallergenic

# Nature's Wealth

# Natural Health Coming From the Past

Did you know that bamboo is the fastest growing plant in the world and can grow up to a yard a day?

Bamboo is a remarkable and versatile natural resource. For centuries, it has found its place in various areas such as cooking, construction, transportation, textile, and medicine in Asia.

Bamboo natural fabric can offer a healthy sleep dynamic for a variety of reasons. Not only are they extremely comfortable, but they also offer lots of health benefits.

## Anti-Allergenic Sleep Experience

If you have allergies, with its fresh texture, bamboo pillow cover is an excellent option for your sleep quality. Thanks to its moisture-absorbing and hypoallergenic properties, you can finally get the sleep you've always been missing.

### - Bamboo fabrics are hypoallergenic.

This means you don't have to worry about allergy-causing bacteria, dust mites and other particles due to its anti-microbial and antibacterial properties.



You sleep comfortably, bamboo fabric takes care of allergies for you.



### - Bamboo fabrics are natural.

Bamboo is a great green resource that is extremely beneficial to the environment.

It's one of the most sustainable resources.

In addition to its soft structure, Bamboo pillow cover will provide a sound sleep thanks to the moisture absorbing feature of the breathing fabric and the comfort it offers.



G

O

↓

G

R

E

E

N



Tested and approved. Skin friendly.

A close-up photograph of an olive branch with several small, green olives. The leaves are dark green and elongated, while the olives are smaller and rounder. The background is slightly blurred, creating a sense of depth.

*Olive Oil*

Let pillow

your  
cares

your  
skin

# Discover a wonderful sleep.

Because Olive Oil pillow cover, which contains the rejuvenating effect of olive oil, ensures that your skin has a balanced moisture during sleep and allows you to discover a wonderful sleep thanks to its perfect softness.

Olive oil infused pillow cover not only offer a natural and sensual feeling when you touch it.

Helps to remove toxins with its antioxidant effect.



Helps to reduce heart problems.



Is a strong antibacterial.



Helps to nourish your skin.



Tested and approved. Skin friendly.



# Youth Code.

C O L L A G E N



# Beautify to you.

Hydra pillow fabrics  
are the perfect  
complement to your  
daily skincare routine.  
Because our skin is  
our first impression.



- Collagen is a very powerful antioxidant that helps to improve the signs of an aging skin and removes the layer of dead cells.



- Collagen is known for its calming, moisturizing and anti-aging effect on all skin types. Its miraculous power helps hydrate the outer layers of your skin in a much softer, smoother and brighter way.



- Collagen is a natural protein that maintains the your skin's strength and restores its elasticity at cellular level.

Tested and approved. Skin friendly.



# Feels So Right





# Feel better with the nature's touch

By transforming fabrics to suit everyday consumer needs, Botanical-based Tencel **cellulose fibers** have become a new definition of sustainability and natural comfort for you.




# Feel the gentle essence of nature.

Built through environmentally responsible processes, domestic TENCEL cellulosic fibers give you a long-lasting, natural, comfortable, and pure sleep environment. With its **Lyocell and Fashionable fibers and softness, TENCEL enhances your sleep quality.**

When you touch TENCEL, you can **feel the refined silky craftsmanship** and see its **natural brightness.**

It gently wraps your body and provides natural comfort for sensitive skin. It helps to absorb moisture more efficiently, so your body feels cool and dry during a comfortable night's sleep.



botanic origin



color retention (TENCEL™ Modal fibers)



moisture management (TENCEL™ Lyocell fibers)



gentle on skin (TENCEL™ Lyocell fibers)



contributes to breathability




unfavorable for bacterial growth (TENCEL™ Lyocell fibers)

A close-up photograph of a woman with brown hair, lying in bed with her eyes closed and a slight smile. She is covered up to her chest with a white sheet. The background is a soft, out-of-focus white, suggesting a clean and comfortable bedroom environment.

# A Comfortable Insect Free Sleep

A healthy and comfortable sleep, away from the insect invasion, awaits you.



- Insects infesting pillows are increasingly becoming a problem around the world. These silent assassins try to feed on human blood at night while an innocent human sleeps peacefully. This causes an ailment that causes itching and scarring.
- Bed bugs like to hide in furniture, beddings, luggages, clothings and are extremely difficult to eradicate once they invade an area.
- Thanks to Healthguard appliqué bed fabric's special formulation creates a warning system against bed bugs, deterring them from infesting the treated item.
- In addition to doing all these, the Healthguard formula also provides a hypoallergenic and hygienic environment.



Tested and approved. Skin friendly.

COOLMAX®

# Dry & Fit Technology


T H E R M A L  
C O M F O R T



# Discover the comfort of cool sleep thanks to the awesome technology of Dry&Fit that is specially produced for you!

We sweat independently from hot and cold weather on a regular basis. Especially during sleep, sweating increases and reduces sleep quality. Consequently, it becomes difficult to sleep in a humid weather.

**Therefore, a mattress made with Dry&Fit Technology improves sleep quality by eliminating the negative effects of sweating.**



## What is Dry&Fit Technology?


The Dry&Fit Technology brand is a family of polyester fibers designed to help you beat the heat. This cooling technology creates a mattress with **permanent moisture absorbing performance**.

Whether it's a hot day or you want to stay cool during sleep, **Dry&Fit Technology helps you stay fresh, dry, and comfortable by evaporating moisture away from your skin and the mattress surface.**

Tested and approved. Skin friendly.

# Albus

FW-317



60x65x20

**SIZE**

## ÖZELLİKLER:

Yatış pozisyonunu ayarladığından mide asidinin yemek borusuna kaçmasını engeller. Reflü rahatsızlığı bulunan kişilere özel olarak tasarlanmıştır. Tam ortopediktir. Sünger bloktan üretilmiştir.

## FEATURES:

Since it adjusts the sleeping position, it prevents stomach acid from escaping into the esophagus. It is specially designed for people with reflux disease. It is fully orthopedic. It is made of sponge block.

## تأثيره

هونلا عضو طببضي منال ارظن نم قدعلا ضمح عنمي مناف ، ءيرملا ىلا بورطلا صاخشال اصي صخ مدمرم من ا رزجلا ضررم نم نواعي نيذلا لمالكاب ماظعل ميوقت من ا قلتك نم ةعونصم يهو جنفسالا



# Carnia

FW-319



SIZE

43X43X9

## ÖZELLİKLER:

Halkanın iç ve dış uçlarına eğim verilerek yere tam oturması sağlanmıştır. Normal oturma minderlerine ek olarak, ön kısmı da boş olmak şartıyla erkek ergonomisine uygun olarak üretilmiştir. Hemoroid ( basur ) , dekubitus yaralarında, kalça ağrılarında ve kosiks dejenerasyonlarında kullanılır.

## FEATURES:

By inclining the inner and outer ends of the ring, it is ensured that it is fully seated on the ground. In addition to the normal seat cushions, it is produced in accordance with male ergonomics, provided that the front part is also empty. It is used in hemorrhoids (hemorrhoids), decubitus wounds, hip pain and coccyx degenerations.

## تأثيره

تيل خادلا فارطال ا قلام! لال خ نم متي ، ققل حليل في ج راخلا او يلع لم الكلاب امتابث نم ذك انتلا دئاسو يلا ففاضال اب. ضرال ا ادج اتنا متي ، في داغلا دغقلا ، فيروكنا لم غلا في يبل افسو يمامالا عزجلا نوكي نا طرشب ممدخ متسا متي. اضيا اغراف (ريس اوبلا) ريس اوبلا يف لكرول مال او اقل تسال جورجو. صن غلا طاطحن او



# Nidus

FW-320



%VISCO&FLEX



SIZE

43x37x7

## ÖZELLİKLER:

Anatomik girintileri sayesinde rahat bir oturma sağlar ve omurliliğin dik durmasına yardımcı olur. Bu sayede uzun süre oturmak zorunda kalan kişilerde bel ve sırt ağrılarının hafiflemesine yardımcı olur. Kalçanızın sizi rahatsız etmeden oturmanızı sağlar. Kaliteli visco malzemeden üretilmiş tam ortopedik bir üründür. Elastik visco malzeme vücut şeklini alır. Basınç ortadan kalkınca yeniden eski haline döner.

## FEATURES:

Thanks to its anatomical recesses, it provides a comfortable sitting and helps the spinal cord to stand upright. In this way, it helps people who have to sit for a long time to relieve lower back and back pain. It allows you to sit without your hips bothering you. It is a fully orthopedic product made of quality visco material. Elastic visco material takes the shape of the body. When the pressure is removed, it returns to its original state.

## تازي مل

تازي مل يورطسي نيدال صاخشالا ايلع تليوط قرتفل سولچلا رهظلا لفسا مالأ فيفخت سولچلاب كل حمسي رهظلاو نل نيكرولا كج عزي نأ نود ماطعلا ميوقتتل لماك جتنم تيلاع وكسييف داوم نم عونصم ذخات قزملل جزللا قدام. قنوجلا طغضلا لاوز دنع مسجلا لكش قيلصالا متلاح ايل دوعي ،

# Coccyx

FW-321



**SIZE**

Flex 44x34x9  
Flex44x34x6  
Visco 44x34x9  
Visco 44x34x6

## ÖZELLİKLER:

Oturma simidini kullandığınızda ilgili bölgedeki ağrı şikayetleri kaybolur veya azalır, basur tedavisinde başarıyı artırır, özellikle kuyruk sokumu ile ilgili enflamasyon semptomlarını yok eder, kalıcı rahatlık sağlar. Kaliteli visco malzemeden üretilmiş tam bir üründür. Elastik visco malzeme vücut şeklini alır.

## FEATURES:

When you use the sitting ring, the complaints of pain in the relevant area disappear or decrease, increase the success in the treatment of hemorrhoids, eliminate the symptoms of inflammation, especially in the coccyx, and provide permanent comfort. It is a complete product made of quality visco material. Elastic visco material takes the shape of the body.

## تأثيرها


، سولجلا ققلا ح مادختسا دنع ،  
عوالكش صرقنت و ايفتخت  
قلصل اناذ ققطنملا يف ملال ا  
يف حاجنلا نم ديزتو ،  
ليزتو ، ريساوبلا جالع  
يف صاخ ، باهتلال اضارعا  
قدمئاد عار رفوتو ، صعصعلا  
داوم نم عونصرم لمك جستنم ن!  
قدم. قنوجلا ايلاع وكسيف  
مسجلا لكش ذخات قنملا جزلل





# Quard Jet

FW-324



**SIZE**

45X43X7

## ÖZELLİKLER:

Ofis sandalyeleri, bilgisayar sandalyeleri, kanepeler, yatar sandalyeler, arabalar, kamyonlar, uçak seyahatleri, sınıflar veya yoga yaparken yere sermek için idealdir.

## FEATURES:

Ideal for office chairs, computer chairs, sofas, recliners, cars, trucks, plane travel, classes or laying on the floor while doing yoga.

## تاليملا

ببتكلملا يساركمل ةيلاشم رتوي بمكلا يساركمو ، دعاقملاو ، لكئارالو ، تانحاشللو ، تارايسلو ، سورنلو ، قرئاطلاب رفسلو ، انشا ضرالو يلوع اقلتسالو و اجويلا قسرامم .



# Hora

FW-326



**SIZE**

50x44x9



## ÖZELLİKLER:

Özel tasarımı sayesinde belinizin bütün boşluklarını tam anlamıyla doldurur. Vücut şeklini alabilen yüksek kalite visco malzemeden üretilmiş olup, tasarımı sayesinde belinizin tüm boşluklarını doldurarak, ekstra destek olur. Uzun süreli oturmalarda rahatlık sağlar.

## FEATURES:

Thanks to its special design, it completely fills all the gaps of your waist. It is made of high quality visco material that can take the shape of the body, and thanks to its design, it fills all the gaps of your waist and provides extra support. It provides comfort during long sittings.

## تاليزملا

ألومي ، صراخلا مديمصرت لضفجب  
ننإ . لكرصخ تاوجف لك أمامت  
ةيلاع وكسيف قدام نم عونصم  
لكش ذخات نأ نكمي قنوجلا  
مديمصرت لضفجبو ، مسجلا  
لكرصخ تاوجف عيمج ألمي  
رفوي . ايفاضل أمعد رفويو  
قليوطلا تاسلجلا ءانثا ةحارلا

# Odina

FW-327



**SIZE**

42x18x10

%100 VISCO

## ÖZELLİKLER:

Baş-boyun-omurga bölgesini dengeler, aynı düzlemde kalmasını sağlar. Boyun fıtığı, boyun düzleşmesi, kireçlenme, duruş ve strese bağlı ağrılara karşı kullanılır. Boyun rahatsızlığı olanların yanı sıra bel ve sırt için de kullanılabilir. Reflü rahatsızlığı olanlar reflü yastığının üzerine ek yapabilir.

## FEATURES:

Balances the head-neck-spine area, keeps it in the same plane. Neck hernia, neck flattening, calcification, It is used against posture and stress-related pain. As well as those with neck problems It can also be used for waist and back. Those with reflux It can add to the reflux pillow.

## تأثيرها

قبقرلاو سارلا ققطنم نزاوي ، يرقفلا دومعلاو ، يوتسمل سفن يف هيقيبى . قنعالا حي طستو قبقرلا قتف سلكنتلاو ملال دض ممدختسا متي ةي عرضول اب طبترملا رتوتلاو نم نوناعي نيذلا كىلوا لىلذكو قبقرلا يف لكاشم رصخلل اضيا ممدختسا نكمي رهظلاو نم نوناعي نيذلا كىلوا عاچترال قداسو ىلا فيضرت نأ نكمي



# Dux

FW-328



**SIZE**

25x19x14

## ÖZELLİKLER:

Dizlerinizin arasına bir yastık yerleştirmek doğru anatomik uyum sağlayarak bacaklar, kalça, diz ve omurgayı konumlandırarak alt bedene destek olur. Bel ağrısı, kas krampları, varisli damarlar ve siyatik hastalıklarından müzdarip kişilerde sorunları azaltmaya yardımcı olur.

## FEATURES:

Placing a pillow between your knees provides correct anatomical alignment and supports the lower body by positioning the legs, hips, knees and spine. It helps reduce problems in people suffering from low back pain, muscle cramps, varicose veins and sciatica.

## تأثيره

ن يجب قداسو عضو رفوي  
ذبحيرشيت ةاذاحم كيتبكر  
ءزجلا معديو ةحيحص  
قيوط نع مسجلا نم يلفسلا  
نيكرولاو نيقياسلا عضو  
يوقفلا دومعلاو نيبتبكرلاو  
لكاشملا ليلقت يف دعاسي  
نوناعي نيذلا صاخشالا يدل  
تاجنشيتو رمظلا لفسا مالآ نم  
قرعو قروالا يلاودو تالضعلا  
رزجلا اسنلا

# Ovis

FW-329



**SIZE**

44x34x9

## ÖZELLİKLER:

Hastalara zarar vermemek ve canlarını yakmamak isteyen refakatçi ve bakıcılar hastayı çevirmeyi sağlayan bu aparatları kullanabiliyor. Hastanın yatakta yatarken bir taraftan bir tarafa dönmesi gerektiğinde hasta çevirme aparatları kullanılıyor.

## FEATURES:

Companions and caregivers who want not to harm or hurt the patients can use these apparatuses that enable the patient to turn. When the patient needs to turn from one side to the other while lying in bed, patient turning apparatuses are used.

## تأثيرها

يهدفون إلى توفير الراحة للمرضى الذين لا يريدون أن يضرهم أو يؤذيهم. يمكن استخدام هذه الأجهزة التي تتيح للمريض التحويل. عندما يحتاج المريض إلى التحويل من جانب إلى الجانب الأخرى أثناء الاستلقاء في السرير، تُستخدم أجهزة تحويل المريض.

# Iter

FW-330



**SIZE**

31x24x10

## ÖZELLİKLER:

Standart şekilde boyun yastığı olarak tasarlanmış ve seyahat sırasında rahatlıkla kullanılabilen bir tercih sağlamaktadır. Bu yastıkları pratik bir şekilde boyun bölgenize yerleştirdikten sonra koltukta oturma sırasında rahatlıkla uyuyabilir ve dinlenebilirsiniz.

## FEATURES:

It is designed as a standard neck pillow and provides a choice that can be used comfortably during travel. After placing these pillows practically on your neck area, you can sleep and rest comfortably while sitting on the sofa.

## تأثيرها

فبقر قئاسوك قءمصم ان! نكءمي ارابخ رفوتو ءيسايق انشا حيرم لكشرب مءءءتسا ءئاسولا مءه عءضو ءعب. رفسلا ، لكشرب ققطنم ءلع اءلمع لكشرب ءءارلا ءونلا لككءمي ءلع سولءلا انشا حيرم ءكءرالا







Feel Well visco orthopedic pillow, with its foam feature, takes the shape of your neck structure when you lie down. The orthopedic pillow, which you can use for neck arthritis, neck hernia, pain caused by fibromyalgia, and neck straightening, can be combined with Feel Well's multi-purpose or cylinder pillow according to your needs.



Supports Neck and Shoulder



Relieves Blood Circulation



Memory Sponge Technology



## Neck Support Pillow

PRODUCT CODE	SIZE	RAW MATERIALS	CASE OPTION
FW-331	60x40x11/8 cm	VISCO	AVAILABLE
FW-332	58x40x11/8 cm	VISCO	AVAILABLE
FW-333	56x36x11/7 cm	VISCO	AVAILABLE
FW-334	51x36x11/7 cm	VISCO	AVAILABLE
FW-335	55x35x13/8 cm	VISCO	AVAILABLE
FW-336	45x31x10/7 cm	VISCO	AVAILABLE



## Cervical Pillow

PRODUCT CODE	SIZE	RAW MATERIALS	CASE OPTION
FW-337	60x34x11 cm	VISCO	AVAILABLE
FW-338	56x34x10.5 cm	VISCO	AVAILABLE



## Neck Support Pillow

DIMENSIONS: 57x42x13/11 cm  
 RAW MATERIAL: VISCO  
 CASE OPTION: AVAILABLE



FW-339



Feel Well, which offers both classic pillow and visco comfort, is preferred with its ideal height for the neck and wide size. The different combinations of this pillow model, with its structure that does not collapse and at the same time feels soft thanks to its dense visco material, appeal to those looking for a high level of comfort.



Supports Neck and Shoulder



Relieves Blood Circulation




Memory Sponge Technology





## Classic Pillow

DIMENSIONS: 60x40x15 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-340

## Contoured Classic Pillow

DIMENSIONS: 60x40x12 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-341

## Oval Classic Pillow


DIMENSIONS: 60x40x12 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE




FW-342

## Jumbo Pillow

DIMENSIONS: 80x40x15 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-343



Feel Well Snoring Reducing Pillow was produced exactly for this need. Feel Well orthopedic pillow, with its structure that fully embraces the neck, prevents possible difficulties in breathing. Suitable for back and side sleeping.



Supports Neck and Shoulder



Relieves Blood Circulation



Memory Sponge Technology



## Snoring Reducer Pillow


DIMENSIONS: 55x40x12/10 cm  
 RAW MATERIAL: VISCO  
 CASE OPTION: AVAILABLE



FW-344

## Apnea Pillow

DIMENSIONS: 55x40x12/10 cm  
 RAW MATERIAL: SPONGE  
 CASE OPTION: AVAILABLE



FW-345

## Butterfly Pillow

DIMENSIONS: 57x45x13 cm  
 RAW MATERIAL: VISCO  
 CASE OPTION: AVAILABLE



FW-346

## Neck Pillow With Arm Support

DIMENSIONS: 65x49x31/12 cm  
 RAW MATERIAL: VISCO  
 CASE OPTION: AVAILABLE



FW-347

## Reflux Pillow

DIMENSIONS: 65x60x24 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-348

## Reflux Pillow

DIMENSIONS: 65x60x24 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-349

## Baby Reflux Pillow

DIMENSIONS: 55x50x10 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-350

## Multipurpose Pillow

DIMENSIONS: 42x18x9 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-351

## Multipurpose Pillow

DIMENSIONS: 42x18x9 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-352

## Cylinder Pillow


DIMENSIONS: 46x14 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-353

## Cylinder Pillow

DIMENSIONS: 46x14 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-354

## Knee Pillow

DIMENSIONS: 25x19x14 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-355

## Patient Turn Pillow

DIMENSIONS: 49x39x13,5 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-356

## Knee Pillow


DIMENSIONS: 25x25x15 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-357

## Varicose Vein Pillow

DIMENSIONS: 61x53x20 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-358

## Oval Bbl Pillow

DIMENSIONS: 41x19x10 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE




FW-359

## Bbl Pillow

DIMENSIONS: 41x19x10 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-360



Feel Well Aroma Therapy Pillow makes it easier for you to fall asleep and provides comfort with its air-channel structure and the fresh scent of mint essence. It was produced for those who prefer classic pillows but want to create a refreshing effect.



Supports Neck and Shoulder



Relieves Blood Circulation



Memory Sponge Technology





## Aroma Therapy Pillows

DIMENSIONS: 60x40x15 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE  
AROMA: LAVANDER



FW-361

## Aroma Therapy Pillows


DIMENSIONS: 60x40x15 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE  
AROMA: MINT



FW-362

## Aroma Therapy Pillows

DIMENSIONS: 60x40x15 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE  
AROMA: OCEAN



FW-363



Our visco orthopedic neck support pillow, supported by cooling gel technology, helps support body temperature. Feel Well Orthopedic Pillow, which can be preferred by those who have sweating problems or want a feeling of coolness in hot weather, is also suitable for those who feel pain in the head and neck area.



Supports Neck and Shoulder



Relieves Blood Circulation



Memory Sponge Technology



## Gel Neck Support Pillow

PRODUCT CODE	SIZE	RAW MATERIALS	CASE OPTION
FW-364	60x40x11/8 cm	VISCO	AVAILABLE
FW-365	58x40x11/8 cm	VISCO	AVAILABLE
FW-366	57x37x11/7 cm	VISCO	AVAILABLE
FW-367	51x37x11/7 cm	VISCO	AVAILABLE
FW-368	55x35x13/8 cm	VISCO	AVAILABLE
FW-369	45x31x10/7 cm	VISCO	AVAILABLE



## Neck Support Pillow

DIMENSIONS: 60x40x15 cm  
 RAW MATERIAL: VISCO  
 CASE OPTION: AVAILABLE



FW-370



It is designed to support your neck and provide comfortable travel on long journeys. With its dense visco texture and different color-fabric options, it protects your neck health during plane, bus and car journeys. Optionally, we also produce it from sponge material.



Supports Neck and Shoulder



Relieves Blood Circulation



Memory Sponge Technology



## Travel Pillow

DIMENSIONS: 31x27x7 cm  
 RAW MATERIAL: VISCO  
 CASE OPTION: AVAILABLE



FW-371

## Travel Pillow

DIMENSIONS: 31x27x10 cm  
 RAW MATERIAL: VISCO  
 CASE OPTION: AVAILABLE



FW-372

## Travel Pillow High

DIMENSIONS: 30x30x12 cm  
 RAW MATERIAL: VISCO  
 CASE OPTION: AVAILABLE



FW-373

## Travel Pillow Large

DIMENSIONS: 31x29x12 cm  
 RAW MATERIAL: VISCO  
 CASE OPTION: AVAILABLE



FW-374



Feel Well's sitting cushions, produced in accordance with the anatomy of the sitting area, also meet different needs. Our cushion, which is used against pain in the sitting area and discomfort caused by sitting on hard floors, has visco and sponge material options.



Relieves Blood Circulation



Memory Sponge Technology



## Sitting Ring

DIMENSIONS: 43x9 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-375

## Sitting Ring

DIMENSIONS: 43x9 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-376

## Sitting Ring

DIMENSIONS: 42x6,5 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-377

## Sitting Ring

DIMENSIONS: 42x6,5 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-378

## Angled Sitting Bench

DIMENSIONS: 43x43x9,5 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-379

## Angled Sitting Bench


DIMENSIONS: 43x43x9,5 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-380

## Wheelchair Cushion


DIMENSIONS: 45x43x7 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-381

## Wheelchair Cushion

DIMENSIONS: 45x43x7 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE




FW-382



## Jelly Square Cushion

DIMENSIONS: 45x43x7 cm  
 RAW MATERIAL: VISCO  
 CASE OPTION: AVAILABLE



FW-383

## Anatomical Cushion

DIMENSIONS: 42x40x9/2 cm  
 RAW MATERIAL: VISCO  
 CASE OPTION: AVAILABLE



FW-384

## Anatomical Cushion


DIMENSIONS: 42x40x9/2 cm  
 RAW MATERIAL: SPONGE  
 CASE OPTION: AVAILABLE



FW-385

## Gel Anatomical Cushion

DIMENSIONS: 42x40x9/2 cm  
 RAW MATERIAL: VISCO  
 CASE OPTION: AVAILABLE



FW-386

## Coccyx Cushion





PRODUCT CODE	SIZE	RAW MATERIALS	CASE OPTION
FW-387	44x36x9,5 cm	VISCO	AVAILABLE
FW-388	44x36x9,5 cm	SPONGE	AVAILABLE
FW-389	44x36x7 cm	VISCO	AVAILABLE
FW-390	44x36x7 cm	SPONGE	AVAILABLE

## Posture Cushion



PRODUCT CODE	SIZE	RAW MATERIALS	CASE OPTION
FW-391	44x38x13/7 cm	VISCO	AVAILABLE
FW-392	44x38x13/7 cm	SPONGE	AVAILABLE





Our orthopedic waist pillow fully supports the waist while sitting and helps to stand upright. Thus, it protects against back, waist and neck pain that may occur in the long term. Our pillow has sponge and visco material options.



Relieves Blood Circulation



Memory Sponge Technology



## Waist Pillow

DIMENSIONS: 35x34x12 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-393

## Waist Pillow

DIMENSIONS: 35x34x12 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-394

## Luxury Waist Pillow

DIMENSIONS: 54x42x9 cm  
RAW MATERIAL: VISCO  
CASE OPTION: AVAILABLE



FW-395

## Luxury Waist Pillow

DIMENSIONS: 54x42x9 cm  
RAW MATERIAL: SPONGE  
CASE OPTION: AVAILABLE



FW-396





FW-399



FW-400

FW-401



FW-402




FW-403



FW-404



FW-404



FW-405







FW-406



FW-407





FW-408



**WALKER**

**30**



FW-409





FW-410



FW-411




FW-412




FW-413





FW-414



FW-415



FW-416



FW-417



WALKER

30



FW-418



FW-419





FW-420



FW-421





FW-422



FW-423



FW-424



FW-425





FW-426



FW-427



FW-428




FW-429



FW-430



FW-431





FW-435



FW-436





FW-437

**OTHER PRODUCTS**

**33**



FW-438





FW-439



FW-440

FW-441



FW-442





**Feel Well**  
MEDICAL



**Feel Well**  
MEDICAL

 [www.feelwellmedical.com](http://www.feelwellmedical.com)